



Clarity AlertMaster Set Up Instructions

1. Open the package, take out all items, and remove all plastic coverings from the items.
2. There should be a total of 5 items in the package:

a. AlertMaster Base Unit



b. Doorbell




c. Bed Shaker



d. Phone Line Cord



e. User Guide



3. Using the cord attached to the Base Unit, plug the system into a wall outlet.
4. To connect to a lamp, first make sure the lamp is turned on, then plug the lamp into the back of the Base Unit.
5. Using the Phone Line Cord, plug one end into the back of the Base Unit, and the other end into the wall. OR plug one end into a TTY/Landline phone and the other into the Base Unit.
6. Plug in the Bed Shaker into the back of the Base Unit and place the shaker under your pillow.
7. Mount the Doorbell outside your door. It should automatically connect after being set up.
8. Set up the clock:
 - First, press the Set Clock button AND the Set hour button at the same time. While holding down the Set Clock button, press the Set hour button until you get to the number you want.

- Then, do the same process to set the minutes, but press the Set minute button instead of the Set hour button.

How to Use:

1. Door light:

- When the Doorbell is pressed, this square, and the sidebars, will light up. The lamp will flash on and off for two minutes. The bed shaker will also vibrate.

2. Phone light:

- When the phone rings, this square will light up, as well as the side bar lights. The lamp will flash on and off for 2 minutes. The bed shaker will also vibrate. To end, answer the phone.



3. Clock light:

- When the alarm clock goes off, this square and the sidebar will light up. The bed shaker will vibrate, and the lamp will flash on



and off for one hour. To end, press the Snooze button.

- When the alarm clock goes off, the bed shaker will vibrate, and the lamp will flash on and off for one hour. To end, press the Snooze button OR the Reset/Page button. *The Snooze is 9 minutes long.

4. Sound light:

- There are 3 modes that the sound indicator can be set to on the back of the Base Unit: All, Main, and Off.
- For ALL: Slide Mic switch to All. When there is a sound that lasts for 12 seconds or longer continuously, the sound square and the sidebars will light up. After another 10

seconds, the bed shaker will vibrate, and the lamp will flash for 2 minutes.

- For MAIN: Slide the Mic switch to Main. When there is a continuous sound for 12 seconds or longer, the sound square and the sidebars will light up. No other indicator will go off.
- For OFF: Slide the Mic switch to Off. No indicator will go off for sound.

****To use the Alarm and Motion indicators, you must purchase additional devices.**

5. To set a Security Timer, use the following buttons:

- Press the Timer on/off button once. This setting turns the lamp on every day at 5pm and off at 11pm.

To confirm this setting, the clock will blink.

- To have the timer set for 6pm, press the Set Timer button and hold for 2 seconds.



- You can do this again to set it for 7pm. After this

setting, it will return to the 5pm setting.

- To turn off the timer for the evening, press the Lamp Control button and hold for 2 seconds. The timer will come on the next evening.

- To turn off the Security Timer completely, press the Timer on/off button once.



1. Snooze button:


- This button is used to temporarily pause an alarm.
The snooze lasts for 9 minutes.

2. Dimmer button:

- This button is used to adjust the brightness of the clock.

3. Reset/Page button:

- This button can be used to stop an alarm clock.
- It can also be pressed to activate any devices that can connect to this system. They will only activate for 2 minutes.



1. Set Clock button:

- Press this button once to set up the time of the clock display.

2. Mic switch:

- This switch is used to control the sound indicator mode.

3. Nightlight Switch:

- This switch is used to turn on and off the nightlight that lights up the ring around the Lamp Control button.